





**Parts Required:** Seat Mount Hanger, Cross Support Tube and Vinyl Strap (configuration specific)  
 -width and depth specific  
 -straight or flared  
 -long or short

**Tools Used:** Scissors  
 3/16" Drill  
 Screw Driver  
 Rubber Mallet

Step	Execution	Figures
<p><b>1</b> A seat mount hanger matrix is available to reference the different options available. See the picture as an example of what the kit includes; two hangers and a cross support tube. Contact your sales representative for configuration verification.</p>		
<p><b>2</b> Disassemble and remove existing legrest and footrest from chair (See legrest removal instructions if required) Make note of the sequence of hardware, this is required when installing the new hangers.</p>		
<p><b>3</b> Remove the front strap from the seat frame by cutting the strap and drilling out the rivets that hold it in place with a 3/16" drill.</p>		
<p><b>4</b> Remove the black plugs at the end of the existing seat frame (this can be done using a screw driver to pry the plugs out)</p>		

- 5** Using a rubber mallet, hammer the hangers into the open tubes in the front of the seat frame.



- 6** Assemble the 2 X 5/16-18 1 1/2" long bolts on each side of the hanger, in the same sequence from step 2.



- 7** Assemble the 2 x 5/16-18 1 1/2" long bolts on each end of the threaded cross support tube to the back of the hangers under the seat.



- 8** Install the vinyl strap supplied (width and color specific) to the front of the chair. See Vinyl Strap Installation instructions if required.



- 9** Attach Swing Away or Elevating Swing Away Footrests.



**For more information, please contact Broda Customer Care at 1-800-668-0637**