






Parts Required: Left & Right Rear Stabilizer

Tools Used: 7/16" Wrench
5/32" Allen Wrench
13mm Wrench

Step	Execution	Figures
1	Rear stabilizers are left and right side specific. Please reference these images for help in determining the correct corresponding side.	 <p style="text-align: right;">Left</p>  <p style="text-align: right;">Right</p>
2	Using a 5/32" allen wrench, remove the MDLC Pedal and control rod. (Located inside the cross bar on the base frame of the chair).	
3	Remove the 6" directional caster using a 7/16" wrench to loosen the bolts and slide the caster out of the caster post. Set caster aside.	

- 4** Slide the stabilizer bar over the caster post tube, then re-insert the caster into the tube frame. Attach the bolt through the hole on the stabilizer frame, the caster post tube and the caster, then tighten to secure.

Note: Set the caster into the directional mode before installing, using an allen wrench or the control rod – see picture below



- 4** Step 4 continued



- 7** Install the Hex Rod through the hole located in the outside wall of the caster post through the casters, leaving approx. the same amount exposed on both sides.



- 8** Install the left and right MDLC Pedals onto the left and right side of the base respectively. Tighten with the 5/32" Allen wrench. Insure that the MDLC pedals are set to the directional mode upon installation. See the corresponding image.



For more information, please contact Broda Customer Care at 1-800-668-0637