

# 587 Seat Height Combinations

Seat Height Post Adjustment Range (With Seat Pan)  (For Strapped Seat Add 1")	Mag Wheel	Mag Wheel	Caster Option 1		Caster Option 2	
	Diameter	Mnt Position	Diameter	Fork Position	Diameter	Fork Position
14", 15", 16", 17"	20"	Mid	6"	High	None	None
15", 16", 17", 18"	20"	High	6"	Mid	8"	High
14", 15", 16", 17"	22"	Low	6"	High	None	None
15", 16", 17", 18"	22"	Mid	6"	Mid	8"	High
16", 17", 18", 19"	22"	High	6"	Low	8"	Mid
15", 16", 17", 18"	24"	Low	6"	Mid	8"	High
16", 17", 18", 19"	24"	Mid	6"	Low	8"	Mid
17", 18", 19", 20"	24"	High	None	None	8"	Low

Note: To achieve the lowest seat height (14" for seat pan, 15" for strapped seat) an extra hole must be drilled in the caster fork (can be done in the factory or in the field) Or a drop seat must be used.