

Where Comfort is Essential



T-8 TRAY

- Slide on dinner tray made of durable ABS plastic
- Great for activities



T8-C CLEAR TRAY

- Slide on clear dinner tray
- Good with Pedal models as it allows residents to see their feet



HALF TRAY

- Durable ABS plastic. Helps provide arm support. Flips down to one side and is fully removable
- Depth adjustable
- Left or right side available
- For push pin height adjustable arms only



T-10 PADDED TRAY COVER

- Provides protection and comfort
- Used with T-8 Tray



MOLDED LATERAL SUPPORTS

- Provides comfortable and durable lateral support and is adjustable along back length and width



S-6 PADDED SEAT BELT

- Pads mid section of resident



S-2 SEAT BELT

- (Side Buckle)



ST-4 PADDED THIGH BELT

- 4" wide
- Supports across thighs instead of mid-section of resident for safety
- Helps prevent resident from sliding



ST-4 NEOPRENE PADDED THIGH BELT

- Machine washable, closed cell water proof foam
- Supports across thighs instead of mid-section of resident for safety
- Helps prevent resident from sliding



SPCV SEAT PAD WITH VINYL COVER PMHBV BACK PAD WITH VINYL COVER

- Durable, easy to clean surface

ESPCV

- Extra Seat Pad Cover- Vinyl

EHBCV

- Extra Back Pad Cover- Vinyl



O2 HOLDER

- Can be mounted on both sides of chair
- O2 Mount required for Midline and Pedal models



N-66 CONTRACTURE SLING

- Offers comfort and support for residents with knee contractures
- Depth adjustable (Model specific)



N-1C HEADREST CUSHION CONTOURED

- (standard foam)



N-1CV HEADREST CUSHION CONTOURED

- (Visco foam)
- Alternative to standard Broda Neckrest for comfort and support



NS NECK SUPPORT

- Half Circle
- Soft Visco foam insert provides comfort around neck
- Cushion rests on resident's shoulders

N-1 NECKREST AND 2 N-7 ROUND BOLSTERS

- Neckrest Cushion Supports head and neck



NI-X-DOUBLE THICK NECKREST

- 5" Thick cushion
- Offers comfort and support and is excellent for residents with kyphosis or hyperextension

Where Comfort is Essential



MAG WHEEL PACKAGE

- 20", 22" or 24" mag wheels (depending on chair model)
- Includes hand rims, anti tippers and wheel locks



N-12 WEDGE BOLSTER

- For various positioning options



HSP PADDING PACKAGE WHEN FACTORY INSTALLED

- Excellent for Huntington's residents
- Includes Short Wing Covers, Extended Shoulder Bolsters, Thigh Belt, Full Extended Footpad and Fixed Footrest
- Pieces also available separately
- Long Wing Covers and Standard Shoulder Bolsters available



HSP PADDING PACKAGE UPGRADE WHEN FACTORY INSTALLED

- Excellent for Huntington's/difficult to seat residents
- Includes standard HSP (left) plus full side panel padding, (3 pc.), and Fleece Legrest Apron and fixed footrest
- Pieces also available separately
- For Elite models
- Long Wing Covers and Standard Shoulder Bolsters available



H-4 IV POLE

- Height adjustable
- Can be mounted on both sides of chair



VENTCART - VENTILATOR CART

- Accommodates ventilator
- Attaches to rear of chair (Elite Models only)
- Ventilator tray available



SAS SIDE ARM SUPPORT

- Covers space between armrest and wing to help keep patient's arms positioned inside the chair



A-5 PADDED ARMREST

- Provides comfort and protection for agitated residents



CUSTOM ARM TROUGH

- Closed cell foam with soft Triptex cover
- Provides support and comfort for resident

CUSTOM LEG TROUGH

- also available



REARSTABIL (Stabilizer)

- Provides additional stability for highly agitated residents

REARSTABLELITE

- (for Elite Models)

FCP-FTMSTABIL FRONT Stabilizer and FCP Lateral Stabilizer also available



F-SW SWING AWAY FOOT RESTS

- (pair) for Pedal models and 587



F-SWELV ELEVATING SWING AWAY ELEVATING FOOTRESTS

- (pair) For Pedal models and 587
- Independently elevating



LEGREST PANEL

- Attaches to swing away footrests to help support resident's legs



FOOTBOARD

- (with Legrest Panel – One piece)
- Attaches to swing away footrests to help support resident's feet and legs



785 XB (EXTENDED BASE)

- Ideal for larger seat depths and heights and for highly agitated residents
- 4.5" forward base extension
- 3" rearward base extension
- 1.5" flare on both sides of front and rear of base (3" total added width)