

Bari 385 Commode/Shower Chair Specifications

Feature	Max	Min	Inc	Comments
Seat Height	23"	21"	2"	Configurable in-field
Seat Width	27"	N/A	N/A	Padded seat width
Effective Seat Width	33"	27"	∞	Distance between armrests
Seat Depth	23"	21"	2"	Configurable in-field
Effective Seat Depth	27"	25"	N/A	Including lower backrest concavity
Seat-Ground Angle	+8°	-14°	∞	+ = Anterior, - = Posterior
Seat Padding	N/A	N/A	N/A	Self-sealing urethane foam
Seat Feature (1)	N/A	N/A	N/A	Open front commode seat hole
Seat Feature (2)	N/A	N/A	N/A	Commode Bed Pan Basket
Backrest Height	26"	N/A	N/A	
Backrest Width	26"	N/A	N/A	
Backrest-Seat Angle	+2°	-65°	∞	+ = Anterior, - = Posterior
Footrest-Seat Distance	20.5"	14"	1.5"	
Footrest Length	10"	N/A	N/A	
Footrest Width	25"	N/A	N/A	
Legrest-Seat Angle	110°	17°	∞	
Armrest-Seat Height	8"	N/A	N/A	
Armrest Length	14"	N/A	N/A	
Armrest Width	3"	N/A	N/A	
Armrest Padding	N/A	N/A	N/A	Molded self-sealing urethane foam
Armrest Feature (1)	N/A	N/A	N/A	Fully removable
Armrest Feature (2)	N/A	N/A	N/A	Adjustable armrests distance (indexed)
Armrest Feature (3)	N/A	N/A	N/A	Armrests tilt with seat
Wheel Configuration	N/A	N/A	N/A	4 locking casters, 2 fixed wheels
Wheel Diameter	8"	5"	N/A	5" locking casters, 8" fixed wheels
Min Turn Radius	N/A	28"	N/A	
Frame Feature (1)	N/A	N/A	N/A	Drainage, no significant water accumulation
Frame Feature (2)	N/A	N/A	N/A	Stainless steel welded construction
Frame Feature (3)	N/A	N/A	N/A	No rusting fasteners or components
Min Ground Clearance	7.75"	2.5"	N/A	Typically 7.75"
Overall Height	52.5"	28.5"	N/A	Typically 48" or 46" in normal seating position
Overall Width	39"	33"	N/A	Typically 33" with armrests retracted
Overall Length	73"	31"	N/A	Typically 46" in normal seating position
Overall Weight	107lbs	N/A	N/A	Without Accessories
Weight Capacity	800lbs	150lbs	N/A	Determined by cylinder lift capacity and testing
ANSI/BIFMA Testing	N/A	N/A	N/A	Meets adapted testing standards