

785 Tilt Recliner

Easily the most advanced chair ever made for long term care, the 785 is a cost effective alternative for hospitals, nursing homes and rehab facilities

List of Features:

1. Up to 90 degrees of back recline (infinitely adjustable)
2. Up to 38 degrees of seat tilt (infinitely adjustable)
3. 16 gauge powder coated tubular steel frame
4. Broda Comfort Tension Seating™
5. Choice of seat or base mounted arms
6. Trendelenburg position capable
7. Leg rest articulates and extends with back recline
8. Leg rest adjustable independently of back recline
9. Leg rest adjusts to any angle between horizontal and vertical
10. Split flip down removable/adjustable footrest
11. Swing away and removable side panels
12. Adjustable Arm Rest Height
13. Removable Wings
14. Removable Padding (shoulder bolsters, neck rest, upper & lower side pads, seat and back pad with removable terry cover, calf pad, sole pad, extra terry covers)
15. 17" seat depth
16. 20.5" standard seat height
17. 18" or 20" seat width
18. 5" total lock front casters, 6" Single directional locking rear casters

Options:

1. Mag wheels (with hand rims, brakes, anti-tip casters)
2. Thigh Belt
3. IV Pole and O2 Holder
4. HSP padding package upgrade
5. Customizing available
6. Depth adjustable armrests
7. Semi recline option
8. Contracture Sling (with semi recline option only)

Intended Users:

Heavy level skilled care facilities
Agitated patients, Ex. Huntington's Disease
Acute hospitals



Elite 785 with HSP padding
and IV Pole