

Broda Pedal Rocker - Model 48R

List of Features:

1. Adjustable Seat Height of 14" – 19"
2. Seat Width 18" or 20"
3. Seat Depth 17"
4. Gentle rocking motion with caregiver lock
5. Up to 18 degrees of infinitely adjustable tilt (33 degrees when combined with rocking function)
6. Adjustable rocking tension to accommodate various residents
7. Broda Comfort Tension Seating™
8. 16 gauge powder coated tubular steel frame
9. Height adjustable armrests
10. Standard removable padding: adjustable neckrest, lower side pads)
11. Non collapsible (anti-theft)
12. 5" Total Lock front casters
13. 5" fixed rear casters
14. Footrest plate brackets

Options:

1. Full padding package: Includes Wings, Shoulder Bolsters, Terry Cloth covered Seat and Back Pad, Extra Terry Covers
2. Thigh Belt
3. Full wrap around vinyl seat and back pad
4. Solid, Clear or Half Tray, Tray padding
5. Mag Wheels
6. IV Pole
7. O2 Holder
8. Swing away or elevating swing away footrests

Intended Users:

Long Term Care Facilities
Residents requiring a light level of care
through to heavy levels of nursing care
Individuals with mild to medium diagnosis
Individuals who are self mobile but lack the
strength or stability



Pedal Rocker with full padding package