

## Broda Acute Care Chair

### List of Features:

- 1 Seat Width 20" or 22"
- 2 Seat Depth 18.5"
- 3 Seat Height 20.5"
- 4 Back Length 32"
- 5 Broda Comfort Tension Seating™
- 6 16 gauge tubular steel frame, powder coated (black) finish
- 7 Up to 27° of infinitely adjustable tilt
- 8 Up to 90° of infinitely adjustable recline
- 9 Removable seat and back for cleaning
- 10 Full padding package- upholstered with minimal and strategically placed stitching to accommodate infection control
- 11 Bolt on wings with removable Shoulder Bolsters
- 12 Swing away, removable, height adjustable arms
- 13 6" Directional locking rear casters
- 14 5" Total locking front casters
- 15 Elevating, articulating and extending legrest with telescopic length adjustment and flip up footrest (16" - 20" range from top of seat)

### Options:

1. Solid or Clear Tray, Tray Padding
2. Mag Wheels (with handrims, wheel locks, anti tippers)
3. IV Pole
4. O2 Holder

**Intended Users:** Acute Care Facilities  
Emergency Wards  
Day Surgery  
Orthopedics  
Long Term Care  
Palliative Care  
Residents/patients requiring medium levels of nursing care



Acute Care Chair  
With optional IV Pole and Tray