

Adult Small Glider Model 100-10 AL

List of Features:

1. 18" seat width
2. 17" seat depth
3. 16" seat height
4. 26" back length
5. Auto Locking Mechanism
6. Broda Comfort Tension Seating™
7. 16 gauge powder coated tubular steel frame
8. Caregiver Lock
9. Soft abrasion-resistant vinyl covered seat and back pad
10. Lumbar Support
11. Rear Transport Wheels
12. Height Adjustable Arms- level to accommodate a tray
13. Comes fully assembled

Options:

1. Clear or Solid Tray
2. Wooden Arms

Intended Users: Long Term Care facilities
Assisted Living facilities
Acute Care
Residents requiring light level of care, or those who are agitated or exhibit a tendency to wander such as residents with Alzheimer's



100-10 AL Adult Small Glider