

## Basic Tilt Recliner- Model 30VT

### List of Features:

1. Up to 45 degrees of infinitely adjustable back recline
2. Up to 25 degrees of infinitely adjustable seat tilt
3. 18" or 20" seat width
4. 17" seat depth
5. Adjustable seat height 19.5" to 21.5"
6. Removable padding (adjustable neck rest, adjustable shoulder bolsters, lower side pads, sole pad)
7. Adjustable armrest height
8. 16 gauge powder coated tubular steel frame
9. Broda Comfort Tension Seating™
10. Legrest articulates with back recline
11. Legrest adjusts independently of back recline
12. Legrest adjusts from any angle between vertical and horizontal
13. Single plate flip up removable/adjustable footrest
14. 3 5" Total lock caster, 1 5" Directional lock rear caster
15. Removable wings

### Options:

1. Full padding package
2. F series split flip up footrest
3. Fixed rear locking casters
4. Solid and/or Clear Trays, Tray Padding
5. Contracture Sling (model would become 30CT)
6. Customizing available

**Intended Users:** General use geriatric chair, superior to traditional geri chairs but less costly than Elite chair

Individuals with mild to medium diagnosis who are generally not self mobile, or who have limited mobility.

Individuals requiring light to medium levels of nursing care



Basic 30VT with full padding, reclined