

## Basic Tilt Chair- Model 80V

### List of Features:

1. Up to 25 degrees of infinitely adjustable seat tilt
2. 16 gauge powder coated tubular steel frame
3. Broda Comfort Tension Seating™
4. Removable padding (adjustable neck rest, adjustable shoulder bolsters, lower side pads, sole pad)
5. 18" or 20" seat width
6. 17" seat depth
7. Adjustable seat height 19.5" to 21.5"
8. Leg rest adjusts independently of seat tilt
9. Leg rest adjusts to any angle between vertical and horizontal
10. Single plate flip down removable/adjustable footrest
11. 3 5" Total lock casters, 1 5" Directional lock rear caster
12. Removable wings

### Options:

1. Wheelchair wheels (with hand rims, brakes, anti-tip bars)
2. Full Padding Package
3. HSP padding package upgrade
4. F series and Split Flip Down footrest
5. Fixed rear locking casters
6. Solid and/or Clear Trays, Tray Padding
7. Contracture sling (model would become 80C)
8. Customizing available

**Intended Users:** Individuals requiring light to medium levels of nursing care

Individuals with mild to medium diagnosis who:

Are generally not self mobile,

Are at risk of slumping or sliding forward,

Who require a fixed angle at the hips

(With wheelchair wheels) individuals who have some self mobility  
but for whom positioning is more important than mobility

Can also be used as Pedal chair by people approximately 6'2" and over.



Basic 80V with full padding