

Commode/Shower Chair – Model CS 385

List of Features:

1. Up to 12 degrees of infinitely adjustable posterior seat tilt
2. Up to 10 degrees of infinitely adjustable anterior seat tilt
3. Up to 67 degrees of infinitely adjustable back recline
4. 22” seat width
5. 18” seat depth
6. 24.5” seat height (to top of padded seat)
7. Broda Comfort Tension Seating™
8. Stainless steel powder coated frame
9. Stainless steel fasteners
10. Swing away and removable arms
11. Removable padded urethane foam commode seat
13. 3 5” Total lock stainless casters
14. 1 5” Directional lock rear stainless caster
15. Single Plate Flip Down removable/adjustable footrest
16. Leg rest adjustable independently of back recline
17. Leg rest adjusts to any angle between horizontal and vertical
18. Commode pan and holder

Options:

1. Adjustable headrest (Padded or Non Padded)
2. Seat belt
5. Side Lateral Supports
6. Extended arms to achieve 25” seat width

Intended Users:

Long Term Care facilities
Acute care facilities
Rehab centers
Individuals who need assistance with shower
and/or commode needs



CS 385 with optional
extended arms