

The Elite Semi-Recliner – Model 35V

List of Features:

1. Up to 45 degrees of infinitely adjustable back recline
2. Up to 25 degrees of infinitely adjustable seat tilt
3. Seat width 18” or 20”
4. Seat depth 17”
5. Back Length 32”
6. 20.5” seat height standard
7. 16 gauge powder coated tubular steel frame
8. Broda Comfort Tension Seating™
9. Removable Padding (shoulder bolsters, upper & lower side pads, neck rest, seat and back pad with removable terry cover, calf pad, sole pad, extra terry covers)
10. Leg rest articulates with back recline
11. Leg rest adjustable independently of back recline
12. Leg rest adjusts to any angle between horizontal and vertical
13. 5” Total lock front casters
14. 6” Single directional lock rear casters
15. Split Flip Down Removable/Adjustable Footrest
16. Swing Away and removable Side Panels
17. Adjustable Armrest Height
18. Removable Wings

Options:

1. Wheelchair wheels (with hand rims, wheel locks, anti-tip casters)
2. HSP Padding Package
3. Contracture Sling (model would become 35C)
4. Solid and/or Clear Trays, Tray Padding
5. Thigh belt
6. Customizing available
7. IV Pole, O2 Holder

Intended Users:

Skilled Long Term Care facilities
Residents requiring medium to heavy levels
of nursing care
Residents who suffer medium to severe diagnosis
needing a repositioning chair with very flexible
positioning capabilities



Elite 35v reclined, with standard padding, wing and side removed