

The Elite Tilt Chair – Model 85V

List of Features:

1. Up to 38 degrees of infinitely adjustable seat tilt
2. 16 gauge powder coated tubular steel frame
3. Broda Comfort Tension Seating™
4. Choice of seat or base mounted arms
5. Removable Padding (shoulder bolsters, neck rest, upper & lower side pads, seat and back pad with removable terry cover, calf pad, sole pad, extra terry covers)
6. Leg rest adjustable independently of seat tilt
7. Leg rest adjusts to any angle between horizontal and vertical
8. 5" Total lock front casters
9. 6" Single directional lock casters
10. Split Flip Up removable/adjustable Footrest
11. Swing Away and removable Side Panels
12. Adjustable Arm Rest Height
13. Removable Wings
14. 20.5 " seat height standard
15. 18" or 20" seat width
16. 17" seat depth

Options:

1. Wheelchair wheels (with hand rims, wheel locks, anti-tip casters)
2. HSP Padding Package
3. Solid and/or Clear Trays, Tray Padding
4. Thigh Belt
5. IV Pole and O2 Holder
6. Contracture Sling (model would become 85C)
7. Depth adjustable armrests
8. Customizing available

Intended Users:

Skilled care nursing homes and other long term care facilities
Residents requiring medium to heavy levels of nursing care
Individuals who suffer medium to severe diagnosis who are not self mobile and need assistance



Elite 85V with base mounted arms