

NEW
Broda Midline Full Recliner

List of Features:

1. Up to 38 degrees of infinitely adjustable seat tilt (varies with seat height)
2. Up to 90 degrees of infinitely adjustable back recline
3. Gas cylinder assisted tilt and recline function
4. 16 gauge powder coated tubular steel frame
5. Broda Comfort Tension Seating™
6. Adjustable armrest height
7. Removable padding (adjustable neckrest, adjustable shoulder bolsters, upper & lower side pads, calf pad, sole pad, seat and back pad with removable terri cover, extra terri covers)
8. 18" or 20" seat width
9. 17" seat depth
10. Adjustable seat height 17" – 22"
11. Legrest elevates with back recline
12. Legrest adjusts independently of back recline to any angle between horizontal and vertical
13. 5" Total Lock front casters
14. 5" Fixed rear casters (or 1 rear 5" Total lock and 1 rear 5" Tracking caster)
15. Single plate flip down removable and length adjustable footrest
16. Removable wings

Options:

1. Cylinder Lock to achieve a semi recline
2. 20", 22" or 24" mag wheels, front or rear mount (with handrims, wheel locks, anti-tippers)
3. HSP Padding Package and Upgraded HSP Padding Package
4. Vinyl Seat and Back Pad Cover
5. F Series Split Flip Down Footrest
6. Solid or Clear Tray, Padded Tray Cover, Solid Half Tray
7. IV Pole, O2 Holder, Ventilator Cart
8. Thigh Belt
9. Lateral Torso Supports
10. Customizing available

- Intended Users:**
- Mid to heavy level skilled care requiring full positioning
 - General resident seating with cylinder lock (for semi recline)
 - Offers more features than the Broda Basic Chair, and is a less expensive alternative with fewer options than the model 785
 - Hospice patients
 - Acute Care



Midline Full Recliner with optional lateral torso supports