

## Broda Pedal Chair™ - Model 48

### List of Features:

1. Adjustable Seat Height of 12.5" – 17.5"
2. Seat Width 18" or 20"
3. Seat Depth 17"
4. Up to 18 degrees of infinitely adjustable seat tilt
5. Broda Comfort Tension Seating™
6. 16 gauge powder coated tubular steel
7. Adjustable Height Armrests
8. Standard Removable Padding: (Adjustable Neck Rest, Lower Side Pads)
9. Care giver operated All Wheel Braking
10. Non collapsible (anti-theft)
11. Rotational 5" front casters
12. Fixed 5" rear casters
13. Footrest plate brackets
14. Care giver or chair user seat tilt adjustment

### Options:

1. Side Wings
2. Shoulder Bolsters
3. Seat and Back Pads
4. Thigh Belt
5. Solid or Clear Trays, Tray Padding
6. Wheelchair wheels (with hand rims, wheel locks, anti tip casters)
7. Seat tilt adjuster mechanism (PAB Bar)
8. User operated brake
9. IV Pole and O2 Holder
10. Customizing available
12. Swing away footrests, or swing away elevating leg rests

### Intended Users:

Long Term Care Facilities  
Residents requiring a light level of care  
through to heavy levels of nursing care  
Individuals with mild to medium diagnosis  
Individuals who are self mobile but lack the  
strength or stability



Model 48 Pedal Chair  
with full padding and clear tray