

Broda 587 Tilt Wheelchair SPECIFICATIONS

Item	Standard	Available Options
Weight Capacity	250 lbs	N / A
Seat Depth	17"	15.5" / 18.5" / 20" w/ strapped seat, 16", 17", 18" 19", 20" w/ seat pan
Seat Height	15" seat pan 16" strapped seat	14" / 16" / 17" / 18" /19"/ 20" seat pan 15" / 17" /18" / 19" / 20" / 21" strapped seat
Seat Width	16" / 18"	20"
Overall Height	42.5" (upright, w/o headrest)	Varies with seat height/wheel options
Overall Width	26.5" / 28.5"	30.5"
Overall Length	35" (w/o Footrests)	Up to 45.5" (w/ Footrests)
Back Height	22"	24", 26", 28", 30"
Degrees of Tilt	Capable of 45° Tilt	
Back Angle	90°	Adjustable to 95°, 100°, 105°
Adjustable Arm Heights (variable, depending on seat height)	5.5" Range	5/8" Increments
Arm Type	Fixed Height, Swing Up Detachable	Height Adjustable, Swing Up, Detachable
Front Casters	6" Swivel	8" Swivel
Rear Wheels	22" Mag Rear Mount	20" or 24" Mag Rear Mount
Frame Color	Black Hammertone	Grey Hammertone
Weight (varies with options)	Approx 66-80 lbs	N / A
Shipping Weight	Approx 75-90 lbs	N / A
Footrest	Swing-Away / Adjustable Length	Swing Away Elevating, Angle Adjustable Footplates
Strap Color	Black	Grey / Blue / Green / Pink