

SPECIFICATIONS



SYNTHESIS



ELITE



MIDLINE



CENTRIC



VANGUARD

| Model number | V4 | 785, 85V | MID | 30VT | 985 |
|---------------------------------|--|--|---|---|---|
| Tilt - Infinitely Adjustable | Up to 40° | Up to 40° | Up to 38° | Up to 25° | 30° Posterior 10° Anterior |
| Recline - Infinitely Adjustable | Up to 90° Trendelenburg capable | Model 785 up to 90° Trendelenburg capable Model 85v-NO recline | Up to 90° Trendelenburg capable | Up to 45° | Up to 45° |
| Frame | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel |
| Standard Padding | Headrest, Shoulder Bolsters, Upper & Lower Side Pads, Calf Pad, Sole Pads, Seat and Back Pad with cover | Headrest, Shoulder Bolsters, Upper & Lower Side Pads, Calf Pad, Sole Pads, Seat and Back Pad with cover | Headrest, Shoulder Bolsters, Upper & Lower Side Pad, Calf Pad, Sole Pads, Seat and Back Pad with cover | Headrest, Shoulder Bolsters, Upper & Lower Side Pads, Calf Pad, Sole Pads, Seat and Back Pad | Headrest, Shoulder Bolsters, Upper & Lower Side Pads, Calf Pad, Seat and Back Pad with cover |
| Upper Lateral Supports | Width and height adjustable frame, removable | Removable | Removable | Removable | Removable |
| Standard Seat Width | 16", 18" or 20" | 16", 18" or 20" | 16", 18" or 20" | 16", 18" or 20" | 26", 28", 30" |
| Standard Seat Depth | 17" or 18.5" | 17" or 18.5" | 17" or 18.5" | 17" or 18.5" | Adjustable 21" or 23" |
| Standard Seat to Floor | 20.5" | 20.5" | Adjustable 17" - 22" | Adjustable 19.5", 20.5", 21.5" | Adjustable 20", 21", 22", 23", 24" |
| Standard Back Length | 32" | 32" | 32" | 32" | Adjustable 32", 34", 36", 38", 40" |
| Custom Seat Widths | 16", 22", 24", 26" | 16", 22", 24", 26" | 16", 22", 24", 26", 28" | 16", 22", 24", 26", 28" | N/A |
| Arm/Side Panels | Height adjustable, swing away and removable | Height adjustable, swing away and removable | Height adjustable | Height adjustable with set screws | Height adjustable |
| Legrest | Length adjustable, Elevates and articulates with back recline. Elevates independently | Length adjustable, Elevates and articulates with back recline. Elevates independently | Length adjustable Elevates independently and with back recline | Length adjustable Elevates independently and with back recline | Length adjustable Elevates independently and with back recline |
| Footrest | Split or single flip down | Split or single flip down | Single flip down | Single flip down | Single flip down |
| Casters | 2 x 5" front total lock 2 x 6" multi-directional | 2 x 5" front total lock 2 x 6" multi-directional | 3 x 5" total lock 1 x 5" directional lock | 3 x 5" total lock 1 x 5" directional lock | 2 x 6" front total lock 2 x 8" multi-directional |
| Weight Capacity | 350 lbs. | 350 lbs. | 350 lbs. | 250 lbs. | 600 lbs. |
| Mag Wheels | Rear or mid mount becomes "550" model | Rear or mid mount becomes "550" model | Rear or front mount becomes "500" model | Rear or front mount becomes "500" model | N/A |
| HCPCS* | N/A | E1161 (785-550SR Only) | N/A | E1031 | N/A |
| ADP Device Code* | N/A | WMA5B-020, WMA5B-010 | N/A | N/A | N/A |

*SPECIFIC CONFIGURATIONS APPLY. CONSULT YOUR REP FOR DETAILS



SASHAY



LATITUDE



COMFORT TILT



TRANQUILLE



ACCESS



REVIVE



REVIVE

| 48 | 48R | 587 | 100-20AL,100-15AL, 100-10AL | LT | CS385 | B385 |
|---|---|--|---|---|---|---|
| Up to 18° | Up to 33° when combined with rocking | Up to 45° | N/A | Up to 28° | Up to 12° Posterior Up to 10° Anterior | Up to 14° Posterior Up to 8° Anterior |
| N/A | N/A | N/A | N/A | N/A | Up to 67° | Up to 65° |
| 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular steel | 16 gauge tubular stainless steel | 16 gauge tubular stainless steel |
| Headrest, Shoulder Bolster, Upper & Lower Side Pads, Seat and Back Pad with cover | Headrest, Shoulder Bolster, Lower Side Pads, Seat and Back Pad with cover | Padding not included (Available as option) | Headrest, Lower Side Pads, Seat and Back Pad (Headrest on Model 100-20AL) | Headrest, seat, back and legrest pad with cover | Padded commode seat | Padded commode seat |
| Removable | Removable | N/A | N/A | N/A | N/A | N/A |
| 16", 18" or 20" | 16", 18" or 20" | 16" or 18" | 18"(100-10AL) 20"(100-15AL,100-20AL) | 19" | 18" or 22" | 27" (33" with arm width adjustment) |
| 17" or 18.5" | 17" or 18.5" | Adjustable 16", 17", 18", 19", 20" | 17" | 18" | 18" | Adjustable 21" or 23" |
| Adjustable 12.5" - 17.5" | Adjustable 14" - 19" | Adjustable 14", 15", 15.5", 16", 16.5", 17", 17.5", 18", 19" | 16"(100-10AL) 18"(100-15AL,100-20AL) | 20" | 24.5"(rolls over most standard toilets) | Adjustable 21" or 23" |
| 32" | 32" | N/A | 26"(100-10AL,100-15AL) 32"(100-20AL) | 32" | 26" | 26" |
| 22", 24", 26", 28" | 22" | 20", 22" | N/A | N/A | Achieve 25" seat width with optional extended arms | N/A |
| Height adjustable | Height adjustable | Removable with Height adjustable option | Height adjustable | Fixed height | Fixed height, swing away & removable | Fixed height & removable |
| N/A Swing away option available | N/A Swing away option available | Swing away option available | N/A | Elevates independently | Length adjustable, Elevated independently of back recline | Length adjustable, Elevated independently of back recline |
| N/A | N/A | Fixed footplate, Angle adjustable footplate available | N/A | N/A | Single flip down | Single flip down |
| 2 x 5" front total lock 2 x 5" rear fixed | 2 x 5" front total lock 2 x 5" rear fixed | 5", 6" or 8" front casters, 20", 22" or 24" rear mag wheels | 2 rear transport wheels | 3 x 5" total lock 1 x 5" directional lock | 3 x 5" total lock stainless 1 x 5" directional stainless | 4 x 7" total lock stainless 2 x 8" middle stainless |
| 350 lbs. | 250 lbs. | 250 lbs. | 250 lbs. | 300 lbs. | 400 lbs. | 800 lbs. |
| Rear or front mount becomes "500" model | Rear or front mount becomes "500" model | Rear mount | N/A | N/A | N/A | N/A |
| E1031 | N/A | E1161 | A9270 | N/A | N/A | N/A |
| WMA5B-005 | N/A | WMA5B-025 | N/A | N/A | N/A | N/A |