



BRODA[®]

Comfort is Essential

COMFORT- TILT
MANUAL WHEELCHAIR

SIMPLIFIED DURABILITY AND VERSATILITY FOR IMPROVED REHAB MOBILITY CARE

The Comfort-Tilt Model 587 provides exceptional durability, positioning and function for complex rehab needs.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 ADJUSTABLE BACK ANGLE
- 3 HEIGHT ADJUSTABLE ARMRESTS
- 4 PADDED REMOVABLE ARMS
- 5 SOLID SEAT PAN
- 6 ADJUSTABLE SEAT TO FLOOR HEIGHT
- 7 ADJUSTABLE SEAT DEPTH

PATIENT BENEFITS



Smooth effortless tilt maintains patient comfort, posture and pressure redistribution.



Optimizes mobility and stability for increased patient safety and function.



Back angle adjustment enhances patient comfort and positioning.



Provides a range of simple adjustments to accommodate a wide variety of needs.

CAREGIVER AND FACILITY BENEFITS

- Posterior tilt enhances pressure redistribution aiding in the maintenance of skin integrity
- Tilt provides ease of use for frequent re-positioning, helping improve patient outcomes
- Innovative design with easy-access fasteners enables quick and simple adjustments
- Proper seating alignment decreases patient fatigue and postural deviations
- Quality construction provides superior durability reducing equipment maintenance and replacement costs

SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Degrees of Tilt	45°
Adjustable Backrest Angle	5°, 10°, 15°, 20°, (13cm, 25cm, 37.5cm, 50cm)
Mag Wheels	20", 22" or 24" (50cm, 55cm, 60cm) Adjustable anti-tippers
Standard Seat Width	16", 18", (40cm, 45cm)
Custom Seat Width	20", 22" (50cm, 55cm)
Standard Seat Depth	16" - 20" (40cm- 50cm)
Standard Seat Height	14" to 19" (35cm-48cm)
Front Casters	5", 6" or 8" with fork (13cm, 15cm, 20cm)
Weight Capacity	250lbs. (114kg)
HCPCS	E1161
ADP Device Code	WMA5B-025



Adjustable seat depth shown

UPHOLSTERY OPTIONS



CUSHION SOLD
SEPERATELY

FRAME OPTIONS



BLACK

COMMON ACCESSORIES

- Swing away footrests
- Elevating swing away legrests
- Angle adjustable footplates
- Heel loops
- Footboard
- Legrest panel
- Seat cushion
- Backrest
- Adjustable headrest
- Lateral torso supports
- Wheel lock extensions
- Aluminum front casters
- Solid or clear tray
- Padded tray cover
- Half tray
- Seat belt
- IV Pole
- O2 Holder



HEALTHCARE DIVISION

info@brodaseating.com

p. 1-844-552-7632

Brodaseating.com

©Broda USA, INC.

8/16