



Smooth and Effortless Caregiver Maneuverability!

Mid-Mount Mag wheels disperse the patient's weight more evenly over the center of the chair base, reducing the front caster load for ease of pushing.

Eliminates lateral drifting for improved control and steering down hallways and turning corners.

Provides excellent tracking on carpeted floors, thresholds and sloped sidewalks/surfaces.

Aids in maneuvering heavier patients, reducing the risk of caregiver injury.

Helps improve overall patient comfort and contentment by offering smooth and efficient maneuverability.

**Mid-Mount Mag wheel package comes with:
14" mag wheels, wheel locks and rear 5" tracking casters
Available factory or non-factory installed
(Shown on Synthesis Tilt Recliner Model V4)**

Industry-Leading Functionality



Rear casters provide anti-tip functionality and are raised 1" from the ground for optimal tracking on uneven surfaces.



Mid-Mount Mags provide excellent maneuverability while tilted or reclined.

Mag wheel position optimizes center of gravity and stability.