



BRODA[®]
Comfort is Essential

REVIVE
BARIATRIC TILT &
RECLINE SHOWER COMMODE

BARIATRIC HYGIENE CARE IN A CLASS ALL ITS OWN

The B385 combines industry leading functionality with BRODA's superior positioning, enhancing patient comfort, care and dignity.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE
- 3 WIDTH ADJUSTABLE, SWING AWAY REMOVABLE ARMS
- 4 4-WAY STRETCH PADDED COMMODE SEAT
- 5 REMOVABLE COMMODE PAN
- 6 INDEPENDENT ELEVATING LENGTH ADJUSTABLE LEGREST
- 7 FLIP DOWN FOOTREST
- 8 MIDDLE 8" WHEEL WITH ONE-TOUCH HANDBRAKE

PATIENT BENEFITS



Posterior tilt improves posture, redistributes pressure and aids in fall prevention.



Back recline accommodates positioning needs for safe and thorough bathing, and access to difficult to reach areas.



Anterior tilt facilitates patient transfers and ease of patient care.



Concave back accommodates gluteal shelf, improving patient support and dignity.

CAREGIVER AND FACILITY BENEFITS

- Enhanced positioning and innovative features facilitates safe transfers
- Provides supportive, stable positioning aiding in exceptional bariatric care
- Commode pan features ergonomic handle and sturdy design for ease of use
- Mid-wheel design enhances maneuverability and prevents lateral drifting, for improved patient transport
- Quality construction provides superior durability reducing equipment maintenance and replacement costs

REVIVE BARIATRIC TILT & RECLINE SHOWER COMMODE B385

COMFORT TENSION SEATING® PROVIDES


- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue



SPECIFICATIONS

Frame	16 Gauge stainless steel frame, with stainless steel components
Backrest	BRODA Comfort Tension Seating®
Seat	Removable, padded commode seat with 4-way stretch cover
Degrees of Tilt	14° Posterior tilt, 8° anterior tilt
Degrees of Recline	65°
Standard Seat Width	27" (68 cm), 33" (83 cm with arm width adjustment)
Overall Width	33" or 39" (83 cm or 98 cm)
Seat Depth	Adjustable 21" or 23" (53 cm or 58 cm)
Seat Height	Adjustable 21" or 23" (53 cm or 58 cm)
Backrest Length	24" (65 cm) stainless steel
Casters	4x5" (12.5) Total lock, 2x8" (20cm) middle wheels, stainless steel
Overall Width	33" or 39" overall width (83 cm or 98 cm)
Weight Capacity	800lbs. (363 kg)

UPHOLSTERY OPTIONS



BLACK
(SEAT ONLY)

FRAME OPTIONS



WHITE

COMMON ACCESSORIES

- Seat belt
- Transport seat
- Molded lateral supports



HEALTHCARE DIVISION

info@brodaseating.com

p. 1-844-552-7632

Brodaseating.com

©Broda USA, INC.

8/16