



BRODA[®]

Comfort is Essential

REVIVE
TILT & RECLINE
SHOWER COMMODE

ENHANCED DIGNITY AND WELL-BEING IN HYGIENE CARE

Superior comfort and versatility, combined with stainless steel frame and components make the CS 385 ideal for shower and commode environments.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE
- 3 HEIGHT ADJUSTABLE HEADREST
- 4 SWING AWAY REMOVABLE ARMS
- 5 4-WAY STRETCH PADDED COMMODE SEAT
- 6 REMOVABLE COMMODE PAN
- 7 INDEPENDENT ELEVATING LENGTH ADJUSTABLE LEGREST
- 8 REMOVABLE FLIP DOWN FOOTREST

PATIENT BENEFITS



Provides positioning for safe and thorough bathing including back recline for washing hair.



Posterior tilt improves posture, redistributes pressure and aids in fall prevention.



Enhances patient comfort and support with full elevating legrest.



Provides anterior tilt and ease of patient care with exceptional maneuverability for positioning over most toilets.

CAREGIVER AND FACILITY BENEFITS

- Enhanced positioning and innovative features facilitates safe transfers
- Provides supportive, stable positioning aiding in exceptional patient care
- Commode pan features ergonomic handle and sturdy design for ease of use
- Quality construction provides superior durability reducing equipment maintenance and replacement costs

COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue




Individualized fit for superior comfort

SPECIFICATIONS

Frame	16 Gauge stainless steel frame, with stainless steel components
Backrest	BRODA Comfort Tension Seating®
Seat	Removable, padded commode seat with 4-way stretch cover
Degrees of Tilt	12° Posterior tilt, 10° Anterior tilt
Degrees of Recline	67°
Standard Seat Width	18" or 22" (46 cm or 56 cm)
Overall Width	24" or 28" (61 cm or 71 cm)
Seat Depth	18" (46 cm)
Seat Height	24.5" (61 cm)
Back Length	26" (66 cm)
Overall Width	24" or 28" (61 cm or 71 cm)
Overall Length	32" (81 cm)
Overall Reclined Length	64" (160 cm)
Front Casters	5" Front total lock, stainless steel (12.5)
Rear Casters	(1) 1x5" Total lock, 1x5" directional lock, stainless steel, (12.5cm)
Weight Capacity	400lbs. (182kg)

UPHOLSTERY OPTIONS



FRAME OPTIONS



WHITE

COMMON ACCESSORIES

- Seat belt
- Transport seat
- Molded lateral supports
- Extended arms (Provide 3" seat width extension)



HEALTHCARE DIVISION

info@brodaseating.com

p. 1-844-552-7632

Brodaseating.com

©Broda USA, INC.

8/16