



**BRODA**<sup>®</sup>  
*Comfort is Essential*

**VANGUARD**  
*TILT RECLINER*

## MORE THAN JUST A WIDER SEAT BARIATRIC CARE RE-DEFINED

Bariatric care is elevated to a new level with the Vanguard Tilt Recliner-Designed specifically for the unique needs of bariatric individuals.



- 1 INFINITELY ADJUSTABLE SEAT TILT
- 2 INFINITELY ADJUSTABLE BACK RECLINE
- 3 REMOVABLE WINGS FOR UPPER LATERAL SUPPORT
- 4 HEIGHT ADJUSTABLE ARMRESTS
- 5 SWING AWAY REMOVABLE ARMS
- 6 STAND-ASSIST GRIPS
- 7 INDEPENDENT ELEVATING AND ARTICULATING LENGTH ADJUSTABLE LEGREST
- 8 FLIP-DOWN FOOTREST
- 9 MULTI-DIRECTIONAL LOCKING REAR CASTERS

### PATIENT BENEFITS



Optional concave back accommodates gluteal shelf, improving patient support and dignity.



Posterior tilt enhances posture, improving digestion, respiratory function and patient endurance.



Manual foot pump-activated seat tilt facilitates re-positioning and transfers, promoting safe patient handling.



Enhanced positioning and innovative features allow for fewer transfers and improved patient comfort.

### CAREGIVER AND FACILITY BENEFITS

- Aids in the maintenance of skin integrity with enhanced pressure redistribution
- Design provides industry leading safety and stability helping to reduce risk of injury
- Tilt and recline provides a change in position for individuals who cannot independently shift their weight
- Quality construction provides superior durability, reducing equipment maintenance and replacement costs

## COMFORT TENSION SEATING® PROVIDES

- Long term seating comfort
- Conforming support for proper pelvic positioning
- Improved air flow, decreasing heat and moisture build up
- Pressure redistribution that aids in maintaining skin integrity
- Accommodation for conditions such as kyphosis and lordosis
- Decreased pressure, improved blood flow and oxygenation of the tissue




CTS Strapping flexes for individualized fit & superior comfort.

## SPECIFICATIONS

Frame	16 Gauge tubular steel frame
Seat & Backrest	BRODA Comfort Tension Seating®
Degrees of Tilt	30° Posterior, 10° Anterior
Degrees of Recline	45°
Standard Seat Width	26", 28", 30"(65cm, 70cm, 75cm)
Standard Seat Depth	21" or 23"(53cm or 58cm)
Standard Seat Height	21", 22", 23", 24"(53cm, 55cm, 58cm, 60cm)
Back Length	32", 34", 36", 38", 40"(80cm, 85cm, 90cm, 100cm)
Front Casters	6" Front total lock (15cm)
Rear Casters	8" Fixed multi-directional locking (20cm)
Weight Capacity	600lbs.(272kg)
Standard Padding	Headrest cushion, upper and lower side pads, shoulder bolsters, seat and back pad with cover, calf pad

## UPHOLSTERY OPTIONS



RENAISSANCE  
BURGUNDY



RENAISSANCE  
CHARCOAL



SILVERTEX  
SAPPHIRE BLUE

## FRAME OPTIONS



GREY  
HAMMER

## COMMON ACCESSORIES

- Tray
- Padded tray cover
- Half tray
- Foot pad
- IV Pole
- O2 Holder
- Thigh belt



HEALTHCARE DIVISION

[info@brodaseating.com](mailto:info@brodaseating.com)

p. 1-844-552-7632

[Brodaseating.com](http://Brodaseating.com)

©Broda USA, INC.

8/16