

MOBILITY CHAIRS



Latitude Pedal Rocker
Model: 48R
Weight Capacity: 250lbs



Sashay Pedal Chair[®]
Model: 48
Weight Capacity: 350lbs

BRODA Comfort Tension Seating[®]
The Key to Pressure Redistribution and Long Term Comfort



LATITUDE PEDAL ROCKER AND SASHAY PEDAL CHAIR[®]
Designed to Encourage Independence and Mobility

Front pivot seat tilt helps prevent falls and allows patient to comfortably self propel, an ideal alternative to a standard wheelchair.

Ergonomic foot on floor position promotes mild exercise.

Adjustable seat to floor height accommodates short and tall patients.

The gentle rocking motion of the Latitude Pedal Rocker provides stimulation and contentment making it excellent for memory care.

Mag wheel option available. (Front of Rear Mount)

POSITIONING CHAIRS



Synthesis Tilt Recliner
Model: V4
Weight Capacity: 350lbs



Midline Tilt Recliner
Model: MID
Weight Capacity: 350lbs

SYNTHESIS TILT RECLINER AND MIDLINE TILT RECLINER
Provides Repositioning for Improved Patient Outcomes

Infinitely adjustable seat tilt and full 90° back recline facilitates proper positioning that improves blood flow, oxygenation and retention of the skin tissue which is vital for maintaining skin integrity.

Full elevating length adjustable legrest with flip-up footplate offers support and improved positioning.

The Synthesis Tilt Recliner offers swing away removable arms for improved transfers and changing of incontinence products.

The Synthesis Tilt Recliner offers an optional WC 19 Compliant Transport Package.

Mag wheel option available.

BARIATRIC POSITIONING



Vanguard Tilt Recliner
Model: 985
Weight Capacity: 600lbs

VANGUARD[®]
TILT RECLINER
Addresses the Unique Needs of Bariatric Patients

Posterior tilt enhances posture, improving digestion, respiratory function and patient endurance.

Anterior Tilt with Stand-Assist Grips facilitate safer transfers.

Optional lower back concavity accommodates patient's gluteal shelf.

Foot pump activated tilt in space provides ease of use for caregivers, reducing the risk of injury.

SHOWER/COMMUNE



Revive Tilt & Recline Shower Commode
Model: CS 385
Weight Capacity: 400lbs

REVIVE TILT & RECLINE SHOWER COMMUNE
Provides Safe Positioning for Bathing

Provides positioning for safe and thorough bathing, including back recline for washing hair.

Provides anterior tilt and ease of patient care with exceptional maneuverability for positioning over most toilets.

Enhances patient comfort, safety and support with full elevating legrest.

Bariatric Model Available up to 800lbs.