

# Comfort Company - Cushions & Backs

Available to order with your Comfort Tilt Chair



www.comfortcompany.com | 800-564-9248

## Cushions:

Curve



HCPCS Coding: E2601 / E 2602

- Mild Anti-thrust shelf
- Medial thigh separator
- High density molded foam

Ascent



HCPCS Coding: E2603 / E2604

- Tapered leg adductors
- Visco ischial insert
- Pelvic, lateral and medial thigh supports

M2



HCPCS Coding: E2607 / E2608

- M2 Zero Elevation
- Tapered leg adductors
- Coccyx relief cutout
- Quadra3D gel pack
- Deep lateral and medial leg contours
- And more...

## Backs:

Acta-Back



Front



Back

- Lightweight
- Provides hip and scapula relief
- 3 layer foam construction
- Laterals and headrests can be attached

Acta-Back Deep



Front



Back

- Lightweight
- Thin profile trunk stabilizer pads
- Provides hip and scapula relief
- Adjustable laterals
- 3 layer foam construction
- Laterals and headrests can be attached



# Comfort Company – Accessories

Available to order with your Comfort Tilt Chair



www.comfortcompany.com | 800-564-9248

## Headrest & Hardware:



- Adjustable: height, depth, angle, side to side
- Secure positioning and attachment
- Moldable to create a custom fit
- Quick-release toggle for easy removal



**Headrest Mounting Hardware**  
HCPCS Coding: E1028

## Lateral Supports:



- Height adjustable
- Width adjustable
- Swing-away or rigid hardware

**Swing-Away Hardware**  
HCPCS Coding: E1028

## Covers:

**STRETCH-AIR** – designed for air transmission and heat reduction. This breathable laminate combines a high-stretch polyester top layer with a breathable vapor barrier

**COMFORT-TEK** – the fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multi-dimensional stretch and helps alleviate pressure on bony prominences

