

Established in 1981, BRODA manufactures high quality tilt in space positioning chairs that provide individuals with the highest level of comfort for the highest quality of care.

By combining BRODA's Comfort Tension Seating<sup>®</sup> with infinitely adjustable tilt and recline, BRODA chairs provide safe, comfortable long term seating that can dramatically reduce the amount of falls that residents face. BRODA chairs work well with "No-Lift" and "Fall Prevention" programs making them the ideal chair to integrate into your safe patient handling practices.

- BRODA offers tilt, recline and legrest adjustments that are operated by gas cylinders. The cylinders help to reduce the amount of lifting force required by caregivers, thus helping to provide safe, frequent repositioning while reducing the risk of injury to the patient and caregiver.
- Both the Synthesis Tilt Recliner and Midline Tilt Recliner provide a full recline position which allows individuals to rest in the chair and allows caregivers to change incontinence products while the individual remains in the chair. This means less transfers!
- The Synthesis Tilt Recliner and other Elite models have removable wings (upper supports) and swing away/removable arms, which allow for safe and effective slide transfers. The footrest flips up accommodating safe stand pivot transfers as well as allowing the use of patient lifts.
- The legrest extends independently with the recline function of the Synthesis Tilt Recliner. The independent elevating legrest helps the individual maintain position as they are being positioned and repositioned throughout the day.
- BRODA chairs work well with most mechanical lifts. When using the deep wedge infinitely adjustable seat tilt, you can accommodate various positioning needs. The combination of tilt and recline accommodates most flexion limitations and postural deviations such as but not limited to, slumping, sliding, falling forward, lateral lean and head drop.
- BRODA chairs provide supportive positioning through a combination of tilt, recline, adjustable legrest angle, wings with shoulder bolsters and height adjustable arms. Supportive positioning decreases postural deviations enhancing patient safety while facilitating safe frequent repositioning throughout the day.
- The Synthesis Transport Chair allows for safe occupant transport in an accessible vehicle with the ability to keep the individual in the same chair for care and transport.