

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____ Ship To Name: _____
 Address: _____ Address: _____
 Phone: _____ Phone: _____
 Contact: _____ Contact: _____
 Email: _____ Tag: _____



Standard Features	
• Comfort Tension Seating®	• 350lb Weight Capacity
• 1 x 5" Rear Directional Locking Caster	• 3 x 5" Total Lock Casters
• Up to 38° Infinitely Adjustable Seat Tilt	• Adjustable Seat to Floor Height
• Up to 90° Infinitely Adjustable Back Recline	• Height Adjustable Armrests
• Tubular Steel, Powder Coated Frame	• Elevating Legrest

Midline Tilt Recliner - Model: MID MSRP \$2,795

Seat and Back Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26" <input type="radio"/> 28"	N/C \$395 \$895	Seat Depth: <input type="radio"/> 15.5" <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"
Seat Height:	<input type="radio"/> 17" <input type="radio"/> 18" <input type="radio"/> 19" <input type="radio"/> 20" <input type="radio"/> 21" <input type="radio"/> 22"	N/C Standard N/C	Back Height: <input type="radio"/> 32" <input type="radio"/> 33.5" <input type="radio"/> 35"
			\$95 Standard N/C \$95 Standard \$95 \$95

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Total Lock Front Casters	Standard	<input type="radio"/> 5" Multi-Directional Locking Rear Casters
<input type="radio"/> 20" Front Mounted Mag Wheels	\$525	<input type="radio"/> 20" Rear Mounted Mag Wheels
<input type="radio"/> 22" Front Mounted Mag Wheels	\$525	<input type="radio"/> 22" Rear Mounted Mag Wheels
<input type="radio"/> 24" Front Mounted Mag Wheels	\$525	<input type="radio"/> Wheel Lock Extension Handles (pair)
		Standard \$525 \$525 \$60

Elevating Footrest Options (Select One)	
<input type="radio"/> Adjustable Elevating Split Footrest (not available in 16" seat width)	Standard
(Choose number of straps) <input type="radio"/> 2 <input type="radio"/> 3 <input type="radio"/> 4 <input type="radio"/> 5	
<input type="radio"/> Adjustable Elevating Single Footrest	N/C
(Choose number of straps) <input type="radio"/> 2 <input type="radio"/> 3 <input type="radio"/> 4 <input type="radio"/> 5	

Armrest Options			
Armrest Options:	<input type="radio"/> Soft Armrest	Standard	Armrest Style:
	<input type="radio"/> Padded Armrests (pair)	\$75	<input type="radio"/> Height Adjustable Arms
	<input type="radio"/> Arm Trough (each)	\$130	<input type="radio"/> Height and Depth Adjustable(pair)
			<input checked="" type="radio"/> Base Mounted Arms
			Standard

Chair Color			
Strap Color:	<input type="radio"/> Grey <input type="radio"/> Blue <input type="radio"/> Black <input type="radio"/> Green	N/C N/C	Padding Color: <input type="radio"/> Burgundy <input type="radio"/> Charcoal <input type="radio"/> Sapphire Blue
Frame Color:	<input type="radio"/> Grey <input type="radio"/> Black	N/C	N/C

Padding - Choose Option 1 OR 2

<input type="radio"/> Option 1 - BRODA Padding				N/C
Standard BRODA Padding:				
Includes: Headrest, Shoulder Bolsters and Wings, Calf & Sole Pad, Upper & Lower Side Pads, Seat & Back Pads with Covers				
Seat & Back Pad Cover	<input type="radio"/> Grey Terry	<input type="radio"/> White Terry		
Options:	<input type="radio"/> Vinyl			N/C
Shoulder Bolster	<input type="radio"/> Standard Short	<input type="radio"/> Rounded Short	<input type="radio"/> ACC Short	N/C
Options:	<input type="radio"/> Standard Long	<input type="radio"/> Rounded Long	<input type="radio"/> ACC Long	\$130
HSP Padding Package:	<input type="radio"/> Includes Wing Cover Long, Thigh Belt & Footbox Long			\$515
	<input type="radio"/> Optional Elevating Fixed Footrest			N/C

 Option 2 - BRODA APP (Additional Positioning Package)

BRODA APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Pads

Seat Width:	<input type="radio"/> 16"	<input type="radio"/> 18"	<input type="radio"/> 20"	Seat Depth:	<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$522
	<input type="radio"/> 22"	<input type="radio"/> 24"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$522
	<input type="radio"/> 26"	<input type="radio"/> 28"			<input type="radio"/> 17"	<input type="radio"/> 19"	<input type="radio"/> 21"	\$522

Common Accessories

Seat Belts:	<input type="radio"/> S-BELT	Seat Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$50	
	<input type="radio"/> S-BELT-A	Seat Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$50	
	<input type="radio"/> S-BELT-PD	Seat Belt Padded			\$72	
	<input type="radio"/> S-BELT-PD-A	Seat Belt Padded Auto Buckle			\$72	
Thigh Belts:	<input type="radio"/> T-BELT	Padded Thigh Belt			\$115	
<i>(not compatible with BRODA APP padding option)</i>	<input type="radio"/> T-BELT-A	Padded Thigh Belt Auto Buckle			\$120	
	<input type="radio"/> T-BELT-NEO	Padded Thigh Belt Neoprene			\$115	
	<input type="radio"/> T-BELT-NEO-A	Padded Thigh Belt Neoprene Auto Buckle			\$120	
Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$295	
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 28" width)	<input type="radio"/> Black Only		\$425	
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap			\$305
	<input type="radio"/> TRAY-HALF	Half Tray (side specific)(includes upper side pad)			\$405	
	<input type="radio"/> TRAY-COVER	Tray Cover (padded)(for TRAY-ABS)			\$150	
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers (select one)(pair)			\$85	
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers			\$150	
Other:	<input type="radio"/> IV-POLE	IV Pole & Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$95	
	<input type="radio"/> O2-HOLDER	O2 Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$160	

Common Padding Upgrades

Headrests:	<input type="radio"/> NR-C	Headrest Contoured Pad	<input type="radio"/> Omit Standard Headrest Less \$60	\$100
	<input type="radio"/> NR-CV	Headrest Visco Foam Contoured Pad		\$100
	<input type="radio"/> NS	Neck Support		\$100
Positioning:	<input type="radio"/> SR	Side Reduction Pads (pair)		\$95
	<input type="radio"/> SAS	Side Arm Supports (pair)		\$125
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports		\$446
	<input type="radio"/> LTS	Lateral Torso Supports (FACTORY INSTALLED RECOMMENDED)		\$550
Footrest Pads:	<input type="radio"/> FB-SH	Footbox - Extended Pad	<input type="radio"/> Omit Calf and Sole Pads Less \$170	\$285
	<input type="radio"/> FB-LONG	Footbox- Full Extended Pad		\$325

Additional Options:
Delivery Instructions or Special Requests:
Total: \$