

BRODA SEATING INC.
When Comfort is Essential
 560 Bingemans Centre Drive, Kitchener, ON Canada N2B 3X9



Phone 1-844-552-7632
 Fax (844)-442-7632
www.brodaseating.com

CDN Ontario MSRP Order Form

Customer Name: _____ Phone #: _____ P.O #: _____

Contact Name: _____ Fax #: _____

Address: _____ Email: _____

Ship to Address: _____ Phone #: _____

Comfort-Tilt MANUAL TILT WHEELCHAIR .. \$2,595 [WMA5B-025](#)

SEAT WIDTH

- 16"..... N/C
- 18"..... N/C
- 20"..... \$195 [WAMA](#)
- 22"..... \$195 [WAMA](#)

SEAT DEPTH

- 16"..... N/C
- 17"..... N/C
- 18"..... N/C
- 19"..... N/C [WAMB](#)
- 20"..... N/C [WAMB](#)

BACKREST ANGLE

- 5°..... N/C
- 10°..... N/C
- 15°..... N/C
- 20°..... N/C

ARM STYLE

- Adjustable Height, Removable..... Included

FRONT RIGGING (select one)

- 75° Swing Away Rigging..... N/C
- 60° Swing Away Rigging..... N/C
- Swing-Away Elevating Legrests..... \$225 [WAMK](#)

FOOTPLATES

- Angle Adjustable, Flip Up Footplates \$100 [WAMJ](#)

FRONT RIGGING HANGER STYLE

- Drop-in (Straight or Flared) (Circle)..... N/C
- Pin and Saddle (Straight or Flared) (Circle)..... N/C

WHEEL LOCK EXTENSIONS

- Wheel Lock Extensions, pair..... included

Standard Features

- Tubular Steel, Powder Coated Frame
- 45 Degrees, Infinitely Adjustable Seat Tilt
- Height Adjustable, Removable Armrests (Arm height adjustment: 9" - 14" from seat)
- Adjustable Anti-Tippers
- Adjustable Seat Depth 16" - 20"
- Adjustable Seat to Floor Height 14" - 19"
- Choice of Wheel, Caster and Fork Sizes
- Push Lock Brakes
- Drop In, Swing Away Foot Rests (4.5" of foot rest length adjustability)
- Angle Adjustable Aluminum Foot Plates
- Padded Side Panels
- Mag Wheel Brake Extensions
- Adjustable Backrest Angle
- 250 lb Weight Capacity

ACCESSORIES

- Slide on Positioning Tray 16" \$ 360
- Slide on Positioning Tray 18"..... \$ 290 [SETND0005](#)
- Slide on Positioning Tray 20", 22"..... \$290 [SETND0010](#)
- Slide on Clear Positioning Tray 18"..... \$350 [SETND0005](#)
- Slide on Clear Positioning Tray 20", 22".... \$ 350 [SETND0010](#)
- Padded Tray Cover..... \$ 170 [SETND0035](#)
- IV Pole..... \$ 180
- O2 Holder..... \$ 200 [WAMX](#)
- Seat Belt 2 point..... \$ 26

FOOTREST OPTIONS

- Heel Loops (pair) \$20
(available only for angle adjustable footplate)
- Legrest Panel (for calf support)..... \$60 [SEFND0010](#)
- Foot Board \$110 [SEFND0015](#)
- Leg Support (one piece)..... \$225 [SEFND0020](#)

BRODA SEATING INC
 WHERE COMFORT IS ESSENTIAL
 560 Binghamans Centre Drive, Kitchener, ON Canada N2B 3X9



Phone 1-844-552-7632
 Fax (844)-442-7632
www.brodaseating.com

CDN Ontario MSRP Order Form

Refer to the Seat Height & Wheel Combinations Matrix and Select:

SEAT TO FLOOR

HEIGHT

- 14"
- 15"
- 15.5"
- 16"
- 16.5"
- 17"
- 17.5"
- 18"
- 19"

REAR WHEELS

- 20" Mag Rear Mount
- 22" Mag Rear Mount
- 24" Mag Rear Mount

FRONT CASTERS

- 5" Front Caster
- 6" Front Caster
- 8" Front Caster

CASTER FORKS

- Small
- Medium

MAG WHEEL POSITION

- A (most forward/maneuverable)
- B
- C
- D* (most rearward/stable)

587V3: SEAT HEIGHT & WHEEL COMBINATIONS								
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB
14"	20	5	SMALL	1	1	TOP	3	2
	20	5	SMALL	2	2	TOP	2	1
15"	22	5	SMALL	1	2	TOP	2	1
	20	6	MEDIUM	2	1	TOP	2	1
15.5"	22	6	MEDIUM	1	1	TOP	2	1
	24	6	SMALL	1	2	TOP	1	1
	24	5	MEDIUM	1	2	TOP	1	1
	24	6	MEDIUM	1	1	TOP	1	1
16"	20	5	SMALL	3	3	TOP	1	1
	20	5	MEDIUM	3	3	TOP	1	1
	20	6	MEDIUM	3	2	TOP	1	1
	22	5	SMALL	2	3	TOP	1	1
	22	6	SMALL	2	2	TOP	1	1
	22	5	MEDIUM	2	3	TOP	1	1
	22	6	MEDIUM	2	2	TOP	1	1
	24	5	SMALL	1	3	TOP	1	1
16.5"	20	6	SMALL	3	3	TOP	1	1
	24	6	SMALL	2	3	BOTTOM	3	1
	24	5	MEDIUM	2	4	BOTTOM	3	1
	24	6	MEDIUM	2	3	BOTTOM	3	1
17"	20	5	MEDIUM	4	5	BOTTOM	3	1
	20	6	MEDIUM	4	4	BOTTOM	3	1
	22	5	MEDIUM	3	5	BOTTOM	3	1
	22	6	MEDIUM	3	4	BOTTOM	3	1
17.5"	20	8	SMALL	4	3	BOTTOM	3	1
	22	8	SMALL	3	3	BOTTOM	3	1
	24	8	SMALL	3	3	BOTTOM	2	1
	24	6	MEDIUM	3	5	BOTTOM	2	1
	24	8	MEDIUM	3	3	BOTTOM	2	1
18"	20	8	MEDIUM	5	4	BOTTOM	2	1
	22	8	MEDIUM	4	4	BOTTOM	2	1
19"	24	8	MEDIUM	4	5	BOTTOM	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.

Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG- 404 587V3 Rev. 1