

PO: _____ **Date:** _____ **Customer #:** _____

Customer Name: _____ Ship To Name: _____
 Address: _____ Address: _____
 Phone: _____ Phone: _____
 Contact: _____ Contact: _____
 Email: _____ Tag: _____



Standard Features	
• Comfort Tension Seating®	• 300lb Weight Capacity
• 16 Gauge Powder Coated Tubular Steel	• 65lbs Overall Weight
• Soft Black Armrests	• 26" Overall Width
• 3 x 5" Total Lock Casters	• Infinitely Adjustable Independent Elevating Legrest
• 1 x 5" Rear Directional Locking Caster	
• Up to 28° Infinitely Adjustable Seat Tilt	

Access Tilt Chair - Model: LT MSRP \$1,300

Seat and Back Options

Seat Width: <input checked="" type="radio"/> 19" Standard	Seat Depth: <input checked="" type="radio"/> 18" Standard
Seat Height: <input checked="" type="radio"/> 20" Standard	Back Height: <input checked="" type="radio"/> 32" Standard

Wheel Options

<input checked="" type="radio"/> 1x 5" Directional Locking Rear Caster & 3x 5" Total Lock Casters	Standard
---	----------

Armrest Options

Armrest Options: <input checked="" type="radio"/> Soft Armrests	Standard
--	----------

Chair Color

Strap Color: <input checked="" type="radio"/> Chestnut Brown Standard	Padding Color: <input checked="" type="radio"/> Meteor Standard
Frame Color: <input checked="" type="radio"/> Chestnut Brown Standard	

Padding

<input type="radio"/> Standard BRODA Padding:	N/C
Includes: Headrest, Legrest Pad, Side Pads, Vinyl Seat & Back pads	
Seat and Back <input type="radio"/> Full Wrap Vinyl Covers (factory installed)	\$106
Pad Cover Options: <input type="radio"/> Seat Pad with V-Notch (for thigh belt)(not available with full wrap vinyl covers)	N/C

Common Accessories

Seat Belts: <input type="radio"/> S-BELT Seat Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$52
<input type="radio"/> S-BELT-PD Seat Belt Padded			\$72
Thigh Belts: (T-Belt only available with v-notch seat pad)			
<input type="radio"/> T-BELT Padded Thigh Belt			\$122
<input type="radio"/> T-BELT-NEO Padded Thigh Belt Neoprene			\$122
Other: <input type="radio"/> IV-POLE IV Pole & Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$100
<input type="radio"/> O2-HOLDER O2 Holder (includes mounting bracket) - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$168

Common Padding Upgrades

Headrests: <input type="radio"/> NR-C Headrest Contoured Pad	<input type="radio"/> Omit Standard	\$106
<input type="radio"/> NR-CV Headrest Visco Foam Contoured Pad	Headrest Less \$64	\$106
Positioning: <input type="radio"/> SR Side Reduction Pads (pair)		\$100

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$ _____