



Customer Care Toll Free: 1-800-668-0637

**Instruction: CYLINDERS BY MODEL**

<b>Document Number:</b>	SWI-07-033
<b>Revision No.:</b>	15
<b>Revision Date:</b>	1-Mar-17
<b>Written By:</b>	Ashley

**Left and Right = When sitting in chair**

Model	Stabulus Part #	Center Mount	Left PAB Cylinder	Right PAB Cylinder	Model	Stabulus Part #	Center Mount	Left PAB Cylinder	Right PAB Cylinder
<b>V4 (Sept 1st 2016 - Current)</b>					<b>385 CS</b>				
Tilt	475856	CYL-800N PAB10			<b>**Switched to Stainless Cylinders March 1st, 2009**</b>				
Recline	10296	CYL-600N PAB17		<b>CYL-600N PAB17 (WC-19)</b>	Tilt	4357		CYLSS-270N PAB 11	CYLSS-270N PAB 11
Legrest	3143LV	CYL-300N PAB5			Recline	4357		CYLSS-270N PAB 11	CYLSS-270N PAB 11
					Legrest	004841		CYLSS-300N PAB 12	
<b>785 (V3) (Sept 16th 2011 - Aug 31st 2016)</b>					<b>550 SR Semi Recline (Option of Full Recline) (July 2016 - Current)</b>				
Tilt - 20" & Greater	475856	CYL-800N PAB10			Tilt - 20" & Greater	475856		CYL-800N PAB10	
Tilt - 16" & 18"	10296	CYL-600N PAB17			Tilt - 16" & 18"	10296		CYL-600N PAB17	
Recline	3143LV		CYL-300N PAB5	<b>CYL-300N PAB5 (WC-19)</b>	Recline	3143LV		CYL-300N PAB5 w/4 x Cyl. Spacer Limiters	
Legrest	3143LV			CYL-300N PAB5	Legrest	3143LV			CYL-300N PAB5
<b>785 (V1 October 2003 - May 2009 &amp; V2 May 2009 - September 2011 ~ Discontinued model)</b>					<b>385 B (Set at 135LBS-lowest weight capacity = Standard)</b>				
Tilt	3186LJ		CYL-270N PAB 4	CYL-270N PAB 4	<b>**Switched to Stainless Cylinders March 1st, 2009**</b>				
Recline	3143LV		CYL-300N PAB5		<b>Note: Option of CYLSS-600N PAB 8 (010296) or CYLSS-800N PAB 9 (010292) - to achieve higher weight - capacity. Please refer to Seat Tilt Lifting Capacity Adjustment Chart in Operating Manual</b>				
Legrest	3143LV		CYL-300N PAB5		Tilt	10297		CYLSS-300N PAB 7	CYLSS-300N PAB 7
<b>35V (Discontinued Model ~ Prior to May, 2009)</b>					Recline	004841		CYLSS-300N PAB 12	CYLSS-300N PAB 12
Tilt	3143LV		CYL-300N PAB5		Legrest	004841		CYLSS-300N PAB 12	
Recline	3143LV		CYL-300N PAB5		<b>Auto Lock Gliders</b>				
Legrest	3143LV			CYL-300N PAB5	Under The Seat	3143LV		CYL-300N PAB 5NL	
<b>85 (V3) (Sept 16, 2011 - Current)</b>					<b>587 V3</b>				
Tilt - 20" & Greater	475856		CYL-800N PAB10		Tilt	691992		CYL-800N PAB 6	
Tilt - 16" & 18"	10296		CYL-600N PAB17		<b>587 V2</b>				
Legrest	3143LV			CYL-300N PAB6	Tilt	691992		CYL-800N PAB 6	CYL-800N PAB 6
<b>85V (V1 &amp; V2 Discontinued)</b>					<b>587 V1</b>				
Tilt	3186LJ		CYL-300N PAB5		Tilt	10297		CYLSS-300N PAB 7	CYLSS-300N PAB 7
Legrest	3143LV		CYL-300N PAB5		<b>Acute Care Chair (Discontinued Model)</b>				
<b>30VT</b>					Tilt	10297		CYLSS-300N PAB 7	CYLSS-300N PAB 7
Tilt	3143LV		CYL-300N PAB5		Recline	3143LV		CYL-300N PAB5	
Recline	3143LV		CYL-300N PAB5		Legrest	3143LV			CYL-300N PAB5
Legrest	3143LV		CYL-300N PAB5		<b>Midline</b>				
<b>80VT (Discontinued Model)</b>					Tilt	3186LJ		CYL-270N PAB 4	CYL-270N PAB 4
Tilt	3143LV		CYL-300N PAB5		Recline	3143LV		CYL-300N PAB5	
Legrest	3143LV		CYL-300N PAB5		Legrest	3143LV		CYL-300N PAB5	
<b>48 Pedal</b>					<b>LT</b>				
Tilt	3186LJ		CYL-270N PAB 4		Tilt	10296		CYL-600N PAB17	
<b>48R Pedal Rocker</b>					Legrest	757398			CYL-200N PAB16
Tilt	21076		CYL-600N PAB 13		<b>Notes:</b>				
<b>985 VG</b>					<b>1. Locking Cylinder includes; cylinder, actuator &amp; lever</b>				
Tilt	MEDACT - H160.1		CYL-HYD-985V2		<b>2. NL = Non Locking (includes; cylinder, actuator &amp; ball bearing)</b>				
Recline	3143LV		2 x CYL-300N PAB5		<b>4. 785, V4, Midline &amp; 30VT Chairs 24" and wider receive one additional Non Locking Recline Cylinder</b>				
Legrest	3143LV			CYL-300N PAB5	<b>5. 48 and 85V Chairs 24" and wider receive one additional Non Locking Tilt Cylinder</b>				
Legrest Length Adjustment Slider			MECHLOK-2	MECHLOK-1					

<b>Cylinder Lengths Measured in Inches</b>			
	<b>Collapsed</b>	<b>Extended</b>	<b>Stroke</b>
<b>PAB 4</b>	10.345	15.4	5.057
<b>PAB 5</b>	12.431	19.73	7.301
<b>PAB 6</b>	12.987	12.75	3.937
<b>PAB 7 (SS)</b>	11.386	17.76	6.378
<b>PAB8 (SS)</b>	11.386	17.76	6.378
<b>PAB 9</b>	11.386	17.76	6.378
<b>PAB 10</b>	11.386	17.76	6.378
<b>PAB 11 (SS)</b>	10.345	15.4	5.075
<b>PAB 12 (SS)</b>	12.431	19.73	7.301
<b>PAB 13</b>	9.775	14.2	4.448
<b>PAB 14</b>	7.516	9.9	2.401
<b>PAB 15</b>	7.516	9.9	2.401
<b>PAB 16</b>	11.2	17.4	6.22
<b>PAB 17</b>	11.386	17.76	6.378