

Enhanced Safety and Comfort

BRODA's HSP Padding Package is designed to provide enhanced patient safety and support while helping to prevent injury.

Our HSP Padding Package provides comfortable, stable, protective padding and positioning for residents with severe agitation and involuntary movement.

Designed to support individuals with:

- Huntington's
- Parkinson's
- ALS (Lou Gehrig Disease)
- Multiple Sclerosis (MS)
- Alzheimer's/Dementia



Synthesis Tilt Recliner Model: V4 with HSP Padding

Factory Installed HSP includes: Long Shoulder Bolsters*

Thigh Belt, Fixed Footrest, Footbox-Long

*Long Shoulder Bolster standard on V4 only, Long Wing Covers standard on other models.



Synthesis Tilt Recliner with up to 40° of infinitely adjustable seat tilt



Synthesis Tilt Recliner with up to 90° of infinitely adjustable back recline

BRODA chairs are the Ideal Solution

BRODA chairs accommodate patients' needs as their conditions change, helping to maintain independence.

Infinitely adjustable positioning cradles and soothes patients, improving posture, reducing agitation and enhancing contentment.

Optimum functionality and durable construction withstand severe agitation and involuntary movements.

BRODA's Comfort Tension Seating[®] paired with gas cylinders, provides energy absorption during episodic movements, enhancing patient safety and comfort.

Lower center of gravity provides industry leading safety and stability aiding in fall prevention.