

Encore

MOBILITY
BRODA[®]

Comfort, Functionality, and Style



Shown with short back, headrest, mag wheels, and lateral support options.



New frame design, including adjustable wings.



Front pivot seat tilt maintains a proper position for self-propulsion.



Pedal Rocker calms and creates contentment.

Patient, Facility, and Caregiver Benefits

- Increases quality of life and engagement with surroundings.
- Allows for easy self-mobility and safer patient maneuverability.
- Reduces risk of injury and falls with supportive stable positioning.



Comfort Tension Seating[®] increases support and sitting tolerance.


The Encore of Comfort, Functionality, and Style.

SPECIFICATIONS


Tilt - Infinitely Adjustable	Rocker: up to 33° (combined with rocking) Chair: up to 25°
Weight Capacity	Rocker: up to 250 lbs. / 158 kg Chair: up to 350 lbs. / 114 kg
Seat Widths Available	Rocker: 16", 18", 20", 22" (40cm, 45cm, 50cm, 55cm) Chair: 16", 18", 20", 22", 24", 26" 28" (40cm, 45cm, 50cm, 55cm, 61cm, 66cm, 71cm)
Seat Depths Available	Rocker: 15.5", 17", 18.5", 20" (39cm, 43cm, 46cm, 50cm) Chair: 17", 18.5", 20" (43cm, 46cm, 50cm)
Back Lengths Available	22", 32", 33.5", 35" (55 cm, 81cm, 85cm, 89cm)
Floor-to-Seat Heights Available	Rocker: Adjustable 14" - 19" (39cm, 43cm, 46cm, 50cm) Chair: Adjustable 12.5" - 17.5"
Wheels - Casters	5" front total lock (12.5cm) 5" rear (12.5cm)
Wheels - Mag Wheels	20" or 22" front mount (50 cm or 55 cm)
Lifetime Warranty Frame	10 years on 16 Gauge tubular steel frame
Pressure-Reducing Seat & Backrest	Integrated with Comfort Tension Seating®
Standard Padding*	Neckrest, shoulder bolsters, lower side pads, seat and back pad with cover

*22" back option excludes neckrest and shoulder bolsters.

UPHOLSTERY OPTIONS



FRAME OPTIONS



Additional accessories available: foot supports, specialty padding packages, trays, and positioning belts