



US MSRP Order Form

Date: _____

Customer Name: _____ Phone #: _____ P.O #: _____

Contact Name: _____ Fax #: _____

Address: _____ Email: _____

Ship to Address: _____ Phone #: _____

HCPCS codes supplied are for reference only and do not ensure coverage or payment. The HCPCS provided are based on our interpretation of the code definitions. HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Comfort-Tilt 587 Manual Tilt Wheelchair \$2,000 E1161

ARMREST

- Fixed Height Armrests N/C
- Height Adjustable, Removable Armrests (factory) ... \$250 E0973

SEAT WIDTH

- 16" N/C
- 18" N/C
- 20" \$195 E2201
- 22" \$195 E2201

SEAT DEPTH

- 16" N/C
- 17" N/C
- 18" N/C
- 19" N/C
- 20" N/C

BACKREST ANGLE

- 5° N/C
- 10° N/C
- 15° N/C
- 20° N/C

FRONT RIGGING (select one)

- 75° Swing Away Rigging..... N/C
- 60° Swing Away Rigging (not available with std. flip up).... N/C
- Swing Away Elevating Legrests \$295 E0990
- F-SW-ELV-IKN Elevating Swing Away
With Angle Adjustable Footplates (pair) \$375

FOOTPLATES (select one)

- Flip up Footplates (on std. swing away rigging only) Standard
- Angle Adjustable, Flip Up Footplates (factory)..... \$185 K0040

FRONT RIGGING HANGER STYLE

- Drop-in (Straight or Flared) (Circle)..... N/C
- Pin and Saddle (Straight or Flared) (Circle)..... N/C
- SM-HANGER-48V-IKN (pair)..... N/C

Standard Features

- Tubular Steel, Powder Coated Frame
- 45 Degrees, Infinitely Adjustable Seat Tilt
- Fixed Height Armrests
- Adjustable Seat Depth 16" - 20"
- Adjustable Seat to Floor Height 14" - 19"
- Choice of Wheel, Caster and Fork Sizes
- Push Lock Brakes
- Drop In, Swing Away Foot Rests
- Composite Foot Plates
- Padded Side Panels
- Adjustable Backrest Angle
- 250 lb Weight Capacity
- Seat cushion, back and headrest not included

ACCESSORIES

- Slide on Positioning Tray 18", 20", 22".... \$290 E0950
- Slide on Positioning Tray 16" \$360 E0950
- Slide on Clear Positioning Tray 18", 20", 22". \$350 E0950
- Padded Tray Cover..... \$170
- IV Pole..... \$180 K0105
- Oxygen Holder..... \$200 E2208
- Seat Belt 2 point..... \$26 E0978
- Heel Loops (pair) \$20 E0951
- Legrest Panel (for calf support)..... \$60 E0995
- Foot Board..... \$110
- Cup Holder \$45
- Accessory Bag (Black) (Each) \$39
- Bag Personalization (2 Initials) Brush Font or \$25
- Block Font

ANTI-TIPPERS

- Rear Anti-Tippers (factory) Required safety feature \$135 E0971

WHEEL LOCKS

- Wheel Lock Extensions, pair..... \$60 E0961



US MSRP Order Form

Refer to the Seat Height & Wheel Combinations Matrix and Select:

SEAT TO FLOOR

HEIGHT

- 14"
- 15"
- 15.5"
- 16"
- 16.5"
- 17"
- 17.5"
- 18"
- 19"

REAR WHEELS

- 20" Mag Rear Mount
- 22" Mag Rear Mount
- 24" Mag Rear Mount

FRONT CASTERS

- 5" Front Caster
- 6" Front Caster
- 8" Front Caster

CASTER FORKS

- Small
- Medium

MAG WHEEL POSITION

- A (most forward/maneuverable)
- B
- C
- D* (most rearward/stable)

*not recommended with 20" wheels

587V3: SEAT HEIGHT & WHEEL COMBINATIONS								
SEAT HEIGHT	MAG WHEEL	CASTER WHEEL	CASTER FORK	MAG PLATE HOLE #	CASTER FORK HOLE #	ANTI-TIPPER MOUNT	ANTI-TIPPER HOLE #	SEAT TAB
14"	20	5	SMALL	1	1	TOP	3	2
15"	20	5	SMALL	2	2	TOP	2	1
	22	5	SMALL	1	2	TOP	2	1
15.5"	20	6	MEDIUM	2	1	TOP	2	1
	22	6	MEDIUM	1	1	TOP	2	1
	24	6	SMALL	1	2	TOP	1	1
	24	5	MEDIUM	1	2	TOP	1	1
16"	24	6	MEDIUM	1	1	TOP	1	1
	20	5	SMALL	3	3	TOP	1	1
	20	5	MEDIUM	3	3	TOP	1	1
	20	6	MEDIUM	3	2	TOP	1	1
	22	5	SMALL	2	3	TOP	1	1
	22	6	SMALL	2	2	TOP	1	1
	22	5	MEDIUM	2	3	TOP	1	1
16.5"	22	6	MEDIUM	2	2	TOP	1	1
	24	5	SMALL	1	3	TOP	1	1
	20	6	SMALL	3	3	TOP	1	1
	24	6	SMALL	2	3	BOTTOM	3	1
17"	24	5	MEDIUM	2	4	BOTTOM	3	1
	24	6	MEDIUM	2	3	BOTTOM	3	1
	20	5	MEDIUM	4	5	BOTTOM	3	1
	20	6	MEDIUM	4	4	BOTTOM	3	1
17.5"	22	5	MEDIUM	3	5	BOTTOM	3	1
	22	6	MEDIUM	3	4	BOTTOM	3	1
	20	8	SMALL	4	3	BOTTOM	3	1
	22	8	SMALL	3	3	BOTTOM	3	1
18"	24	8	SMALL	3	3	BOTTOM	2	1
	24	6	MEDIUM	3	5	BOTTOM	2	1
	24	8	MEDIUM	3	3	BOTTOM	2	1
	20	8	MEDIUM	5	4	BOTTOM	2	1
19"	22	8	MEDIUM	4	4	BOTTOM	2	1
	24	8	MEDIUM	4	5	BOTTOM	1	1

Hole numbering begins at the top of the plate, fork and anti-tipper. Hole number 1 is the uppermost hole.

Do not use the rear most vertical row/column of holes on the mag plate for 20" mag wheels

Form08CONFIG- 404 587V3 Rev. 1

BRODA SEATING INC.

When Comfort is Essential

560 Bingham Centre Drive, Kitchener, ON Canada N2B 3X9



Phone: 1-844-552-7632

Fax (844)-442-7632

www.brodaseating.com

US MSRP Order Form

Comfort Company Items will Ship Separate from Comfort-Tilt



Seat cushion, back, lateral supports and headrest sold with chair only

SEAT CUSHION		16"-20"W 22" W	SEAT CUSHION WIDTH		
<input type="checkbox"/>	Comfort Company Curve®.....	\$88 \$185	E2601	<input type="checkbox"/> 16"	N/C
<input type="checkbox"/>	Comfort Company Ascent	\$245 \$306	E2603	<input type="checkbox"/> 18"	N/C
<input type="checkbox"/>	Comfort Company M2 Zero Elevation	\$390 \$418	E2607	<input type="checkbox"/> 20"	N/C
				<input type="checkbox"/> 22"	N/C
SEAT CUSHION DEPTH			SEAT CUSHION COVER		
<input type="checkbox"/>	16"	N/C	<input type="checkbox"/>	STRETCH-AIR	N/C
<input type="checkbox"/>	18"	N/C	<input type="checkbox"/>	COMFORT-TEK	N/C
<input type="checkbox"/>	20"	N/C			

BACK		16"-20"W 22" W	BACK WIDTH		
<input type="checkbox"/>	Comfort Company Acta-Back® 16" Tall	\$507 \$596	E2613	<input type="checkbox"/> Fits 15" - 17" W/C	N/C
<input type="checkbox"/>	Comfort Company Acta-Back® 18" Tall	\$507 \$596	E2613	<input type="checkbox"/> Fits 17" - 19" W/C	N/C
<input type="checkbox"/>	Comfort Company Acta-Back® 20" Tall	\$576 \$653	E2613	<input type="checkbox"/> Fits 19" - 21" W/C	N/C
				<input type="checkbox"/> Fits 22" W/C	N/C
DEEP CONTOUR BACK		16"-20"W 22" W	BACK CUSHION COVER		
<input type="checkbox"/>	Comfort Company Acta-Back® Deep 16" Tall.....	\$601 \$677	E2615	<input type="checkbox"/> STRETCH-AIR	N/C
<input type="checkbox"/>	Comfort Company Acta-Back® Deep 18" Tall.....	\$601 \$677	E2615	(STRETCH-AIR not available with 22" Tall Back)	
<input type="checkbox"/>	Comfort Company Acta-Back® Deep 20" Tall.....	\$678 \$754	E2615	<input type="checkbox"/> COMFORT-TEK	N/C
<input type="checkbox"/>	Comfort Company Acta-Back® Deep 22" Tall.....	\$678 \$754	E2615		

LATERAL SUPPORTS		LATERAL SUPPORT PAD COVER			
<input type="checkbox"/>	Comfort Company Single Lateral Pad 6"W x 4"T	\$79	E0956	<input type="checkbox"/> STRETCH-AIR	N/C
<input type="checkbox"/>	Comfort Company Single Lateral Pad 7"W x 5"T	\$79	E0956	<input type="checkbox"/> COMFORT-TEK	N/C
<input type="checkbox"/>	Comfort Company Single Rigid Hardware -Left (Hardware Only)	\$125			
<input type="checkbox"/>	Comfort Company Single Rigid Hardware -Right (Hardware Only)	\$125			
<input type="checkbox"/>	Comfort Company Single Swing Away Hardware Left (Hardware Only)	\$263	E1028		
<input type="checkbox"/>	Comfort Company Single Swing Away Hardware Right (Hardware Only)	\$263	E1028		

HEADREST		HEADREST PAD COVER			
<input type="checkbox"/>	Comfort Company BodiLink Headrest Pad, 8"W X 4"L	\$ 230	E0955	<input type="checkbox"/> STRETCH-AIR	N/C
<input type="checkbox"/>	Comfort Company BodiLink Headrest Pad, 10"W X 5"L	\$ 230	E0955	<input type="checkbox"/> COMFORT-TEK	N/C
<input type="checkbox"/>	Comfort Company Removable Headrest Mounting Hardware (Headrest pad not included)	\$ 221	E1028		