

PO: _____

Date: _____

Customer #: _____

Customer Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Email: _____

Ship To Name: _____
 Address: _____
 Phone: _____
 Contact: _____
 Tag: _____



Standard Features	
• Up to 14° Infinitely Adjustable Seat Tilt (Posterior)	• 800lb Weight Capacity
• Up to 8° Infinitely Adjustable Seat Tilt (Anterior)	• Contoured Backrest
• Up to 65° Infinitely Adjustable Back Recline	• Swing Away Removable Armrests
• Stainless Steel, Powder Coated Frame	• Comfort Tension Seating®
• 4 x 6" Total Lock Stainless Steel Front & Rear Casters	
• 2 x 8" Locking Middle Wheels With One-Touch Emergency Brake	

Revive Bariatric Tilt & Recline Shower Commode - Model: B385 MSRP \$5,995

Seat and Back Options			
Seat Width:	<input checked="" type="radio"/> 27" N/C	Seat Depth:	<input type="radio"/> 21" N/C
	(33" with width adjustable armrests)		<input type="radio"/> 23" N/C
Seat Height:	<input type="radio"/> 21" N/C	Back Height:	<input checked="" type="radio"/> 26" Standard
	<input type="radio"/> 23" N/C		

Wheel Options	
<input checked="" type="radio"/> 6" Total Lock Stainless Steel Front & Rear Casters & 8" Locking Middle Wheels with One Touch Emergency Brake	Standard

Seat Options (Select one)	Bed Pan & Basket Options
<input type="radio"/> Seat Commode Padded	Standard
<input type="radio"/> Seat Transport Padded	N/C
	<input checked="" type="radio"/> Bed Pan with Basket
	Standard

Armrest Options	
Armrest Options:	<input checked="" type="radio"/> Soft Armrest Standard
Armtop Style:	<input checked="" type="radio"/> Seat Width Extending (extends width to 33") Standard

Headrest and Lateral Support Options	
Headrest:	<input type="radio"/> Molded 2" Offset \$150

Strap Color	
Strap Color:	<input type="radio"/> Grey <input type="radio"/> Blue <input type="radio"/> Green N/C

Preset Factory Installed Weight Setting (Select one)			
<input type="radio"/> 135 lb	<input type="radio"/> 270 lb	<input type="radio"/> 360 lb	<input type="radio"/> 405 lb N/C
<input type="radio"/> 540 lb	<input type="radio"/> 630 lb	<input type="radio"/> 720 lb	N/C

Common Accessories	
Lateral Supports:	<input type="radio"/> Molded Laterals (pair) \$75
Seat Belt:	<input type="radio"/> Seat Belt (Side Buckle) (Large Size) \$50

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$