

PO: _____

Date: _____

Customer #: _____

Customer Name: _____

Ship To Name: _____

Address: _____

Address: _____

Phone: _____

Phone: _____

Contact: _____

Contact: _____

Email: _____

Tag: _____



Standard Features	
• Comfort Tension Seating®	• 350lb Weight Capacity
• Up to 40° Infinitely Adjustable Seat Tilt	• Trendelenburg Capable
• Up to 90° Infinitely Adjustable Back Recline	• Elevating and Articulating Legrest
• 2 x 6" Rear Multi-Directional Locking Casters	• 2 x 5" Total Lock Front Casters
• Height Adjustable, Swing-Away Removable Arms	
• Height and Width Adjustable Upper Laterals (Wings)	

Synthesis Tilt Recliner - Model: V4 MSRP \$3,674

Seat and Back Options			
Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" <input type="radio"/> 22" <input type="radio"/> 24" <input type="radio"/> 26"	N/C \$418 \$946	Seat Depth: <input type="radio"/> 17" <input type="radio"/> 18.5" <input type="radio"/> 20"
			Standard N/C \$100
Seat Height:	<input type="radio"/> 19.5" <input type="radio"/> 20.5" <input type="radio"/> 22" <input type="radio"/> 23.5"	\$180 Standard \$180	Back Height: <input type="radio"/> 32" <input type="radio"/> 33.5" <input type="radio"/> 35"
			Standard \$100 \$100

Front Options	Wheel Options	Rear Options
<input type="radio"/> 5" Total Lock Front Casters	Standard	<input type="radio"/> 6" Multi-Directional Locking Rear Casters
<input type="radio"/> 5" Aluminum Non-Locking Front Casters	\$190	<input type="radio"/> 14" Rear Mags <input type="radio"/> 14" Mid Mags
<i>(available with rear mags only)</i>		<input type="radio"/> 20" Rear Mags
<input type="radio"/> 5" Front Locking Twin Wheel Casters	\$190	<input type="radio"/> Wheel Lock Extension Handles (pair)
		Standard \$332 \$332 \$64

Elevating Footrest Options (Select One)	
<input type="radio"/> Adjustable Elevating Split Footrest (not available in 16" seat width)	Standard
Choose Number of Straps (footrest length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")	
<input type="radio"/> Adjustable Elevating Single Footrest	N/C
Choose Number of Straps (footrest length range) <input type="radio"/> 2 (12"-19") <input type="radio"/> 3 (14"-21") <input type="radio"/> 4 (16"-23") <input type="radio"/> 5 (18"-25")	
<input type="radio"/> ABS Elevating Leg & Footrest (Includes calf and sole pads) (16" - 24" Widths)	\$344

Armrest Options			
Armrest Options:	<input type="radio"/> Soft Armrest	Standard	Armrest Style: <input type="radio"/> Height Adjustable Arms
	<input type="radio"/> Long Armrest	\$42	<input type="radio"/> Height and Depth Adjustable (pair)
	<input type="radio"/> Padded Armrests (pair)(black)	\$80	<input checked="" type="radio"/> Seat Mounted Arms
	<input type="radio"/> Arm Trough (each)	\$136	Mount:
			Standard

Base Frame Finish Color (Select one)			
<input type="radio"/> Carmine Red	<input type="radio"/> Ocean Blue	<input type="radio"/> Midnight Black	N/C
<input type="radio"/> Satin Nickel	<input type="radio"/> Tungsten Sand	<input type="radio"/> Bronze	N/C

Upholstery Color Options (Select one)

- * Standard padding configuration: Headrest, exterior of shoulder bolsters, upper and lower side pads & sole pads Come in Accent Color. Seat, back and calf pads, and interior of shoulder bolsters come standard in Savoy Black Vinyl.
- ** Indigo, Espresso & Graphite are available on seat, back & calf pad & most accessories. See below under padding for selection.
- *** In honor of our Veterans Service, Camo options are now available.



Standard Padding*: Shown with Geyser Upholstery & Tungsten Sand Frame Color



Accent Padding**: Shown with Indigo Upholstery & Black Frame Color *Chair Shown with ABS footrest

Upholstery Accent Color* (Select One)

Vinyl:

- Indigo**
- Interlace Espresso**
- Graphite**
- Black
- Smoke House
- Chocolate

Polyurethane:

- Portobello
- Geyser
- Fossil
- Forge

Camo Nylon:* (with performance finish)**

- Digital Camo Green
 - Retro Camo Jungle
- \$130

Padding - Choose Option 1 OR 2

Option 1 - BRODA Padding

N/C

Standard BRODA Padding:

Includes: **In Accent Color:** Headrest, Exterior Shoulder Bolsters, Upper & Lower Side Pads & Sole Pads*.

In Black Vinyl: Seat & Back Pads with Covers, Calf Pad & Interior Should Bolsters

Seat, Back & Calf	<input type="radio"/> Black Vinyl	Standard
Pad Cover Options:	<input type="radio"/> Accent Color*	N/C
	*** <i>(For Indigo, Espresso or Graphite only)</i>	
	<input type="radio"/> Full vinyl Wrap Seat & Back Pads	\$100
	<i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	
	<input type="radio"/> Calf Pad Color	N/C
	<i>(Select Black Vinyl OR Accent Color Vinyl Above)</i>	
	<input type="radio"/> Grey Terry <input type="radio"/> White Terry	N/C
HSP Padding Package:	<input type="radio"/> Includes Shoulder Bolsters Long, Thigh Belt & Footbox Long	\$544
	T-Belt Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	Footbox Long Color: <input type="radio"/> Black Vinyl <input type="radio"/> Accent Color Vinyl	N/C
	<input type="radio"/> Optional Elevating Fixed Footrest	N/C

Option 2 - BRODA APP (Additional Positioning Package)

BRODA APP Padding Package Includes: Seat Cushion, Back and Head Support, Foot & Calf Support, Upper & Lower Side Pads

Seat Width:	<input type="radio"/> 16" <input type="radio"/> 18" <input type="radio"/> 20" (E2601, E2611)	Seat Depth:	<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522
	<input type="radio"/> 22" <input type="radio"/> 24"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522
	<input type="radio"/> 26"		<input type="radio"/> 17" <input type="radio"/> 19" <input type="radio"/> 21"	\$522

Acute Care Padding Package:

Includes: Full Vinyl Covers, ABS Footrest with Calf and Sole Pads *(Factory installed)(skip footrest and padding package section above when selecting this option)(16"-24" widths)* \$476

Common Accessories

Seat Belts:	<input type="radio"/> S-BELT	Seat Belt	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$52
	<input type="radio"/> S-BELT-A	Seat Belt Auto Buckle	<input type="radio"/> Center Buckle	<input type="radio"/> Side Buckle	\$52
	<input type="radio"/> S-BELT-PD	Seat Belt Padded			\$76
	<input type="radio"/> S-BELT-PD-A	Seat Belt Padded Auto Buckle			\$76
Thigh Belts: <i>(not compatible with BRODA APP padding option)</i>	<input type="radio"/> T-BELT	Padded Thigh Belt			\$122
	<input type="radio"/> T-BELT-A	Padded Thigh Belt Auto Buckle			\$126
	<input type="radio"/> T-BELT-NEO	Padded Thigh Belt Neoprene			\$122
	<input type="radio"/> T-BELT-NEO-A	Padded Thigh Belt Neoprene Auto Buckle			\$126

Common Accessories (Continued)

Trays:	<input type="radio"/> TRAY-ABS	Tray ABS with Strap	<input type="radio"/> Sewn Tray Strap	<input type="radio"/> Grey <input type="radio"/> Black	\$312
	<input type="radio"/> TRAY-ABS-OVR	Tray ABS with Strap (24" - 26" width)		<input type="radio"/> Black Only	\$448
	<input type="radio"/> TRAY-CLEAR	Tray Clear with Strap	<input type="radio"/> Sewn Tray Strap		\$322
	<input type="radio"/> TRAY-HALF	Half Tray <i>(side specific)(includes upper side pad)</i>			\$428
	<input type="radio"/> TRAY-COVER	Tray Cover <i>(padded)(for TRAY-ABS)</i>			\$158
Stabilizers:	<input type="radio"/> STABIL-FR	<input type="radio"/> Front or <input type="radio"/> Rear Mounted Stabilizers <i>(select one)(pair)</i>			\$90
	<input type="radio"/> STABIL-LAT	Lateral Frame Stabilizers			\$158
Other:	<input type="radio"/> IV-POLE	IV Pole & Holder <i>(includes mounting bracket)</i> - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$100
	<input type="radio"/> O2-HOLDER	O2 Holder <i>(includes mounting bracket)</i> - <input type="radio"/> Front or <input type="radio"/> Rear Mount			\$168
	<input type="radio"/> ACCESSORY-BAG	Accessory Bag (Black) Each			\$39
	<input type="radio"/> BAG PERSONALIZING (2 Initials)	<input type="radio"/> Brush Font OR <input type="radio"/> Block Font			\$25
	<input type="radio"/> CUP-HOLDER	Cup Holder			\$45

Common Padding Upgrades

Headrests:	<input type="radio"/> NR-CV	Headrest Visco Foam Contoured Pad	<input type="radio"/> Omit Standard	\$106
	<input type="radio"/> NS	Neck Support	Headrest Less \$64	\$106
Positioning:	<input type="radio"/> SR	Side Reduction Pads <i>(pair)</i>		\$100
	<input type="radio"/> CC-LTS	Comfort Company Lateral Supports		\$472
	<input type="radio"/> LTS	Lateral Torso Supports <i>(FACTORY INSTALLED RECOMMENDED)</i>		\$580
	<input type="radio"/> SB-LONG-V4	Shoulder Bolster Long <i>(pair)</i> - Omits standard short		\$286
Footrest Pads:	<input type="radio"/> FB-LONG	Footbox- Full Extended Pad	<input type="radio"/> Omit Calf and Sole	\$344

Additional Options: _____

Delivery Instructions or Special Requests: _____

Total: \$
